E n g l i s h	2
F r a n ç a i s	6
D e u t s c h	10
E s p a ñ o l	14
I t a l i a n o	18
P o r t u g u ê s	22

Original Harvey
and Harvey DLX
Double and
Triple Head
Stethoscopes

Designed to enable the medical professional to obtain superior and accurate auscultation.

The Harvey and Harvey DLX Double Head and Triple Head Stethoscopes were designed after extensive research and testing. The critical design of the tubing, proper weight and sizing of the chestpiece and rotatable binaurals, combine to make the Harvey Stethoscopes exceptional instruments for auscultation. The final testing, of course, was on the patient.

The Harvey Stethoscopes will give you the ability to identify the lowest frequencies that the human ear can detect, such as a faint diastolic gallop or rumble, as well as the higher frequencies of the faint early blowing aortic diastolic murmur.

Furthermore, with over five decades of experience using stethoscopes, Dr. W. Proctor Harvey has designed his stethoscopes to meet the comfort and portability needs which physicians demand.

Double Head

Triple Head

Chestpieces

Both the double head and triple head offer a flat diaphragm and bell. The triple head also utilizes a corrugated diaphragm for additional auscultatory capabilities. Selection is accomplished by holding the chestpiece stem in one hand and rotating the chestpiece until a click is heard or felt.

Flat Diaphragm

“The workhorse.” Excellent for all sounds and murmurs. Best for higher frequencies (such as faint diastolic murmur of aortic regurgitation). Also ideal to identify splitting of sounds, systolic clicks and ejection sounds. 4.76 cm diameter.

Corrugated Diaphragm

Amplifying quality. Excellent for an overview of heart sounds. Especially good for low-frequency gallop sounds and murmurs. 4.76 cm diameter.

Bell



Provides excellent detection and clarity for low and medium-frequency heart sounds and murmurs (such as a faint gallop or diastolic rumble). 2.79 cm diameter.

Headset Adjustment


Your new Harvey Stethoscopes are designed to give you the maximum ability to adjust for your personal comfort and acoustic fit.


- A.** Select the correct size and texture eartip. For maximum acoustic seal, two sizes of “soft” and “hard” eartips are provided. Choose the size which best seals off your outer ear canal and feels the most comfortable. As a rule the larger size is best for the majority.

- B.** To adjust the angle of the binaurals, hold the spring in one hand and rotate each binaural to the correct angle.


- C.** Make sure that the binaural tubes are pointed slightly forward going into the ear canal.

D. To reduce spring tension, grasp each eartube at its bend and gently pull outwards.

E. To increase spring tension, squeeze the eartubes together.


Care and Maintenance


- It is important to clean and disinfect your stethoscope to prevent possible cross-contamination.
- Eartubes and tubing can be wiped down with alcohol or mild soapy water.
- Eartips may be removed for thorough cleaning with alcohol or mild soapy water.
- Diaphragm replacement: Under normal use you should rarely need to remove the rim and diaphragm for cleaning purposes. The exterior of the diaphragm can be cleaned with alcohol or mild soapy water. However, if you do need to remove diaphragm, follow the instructions below:

CAUTION: The original Harvey stethoscopes contain rubber latex which may cause allergic reactions.


A. Hold the bell and turret section in one hand, and exert pressure on the back of the diaphragm. Grip the threaded ring of the diaphragm with your other hand and turn the ring counterclockwise to loosen.

B. Separate the threaded ring from the turret section.

C. Simply press the diaphragm disc out of the metal ring and replace with new disc. Thread the ring back on to chestpiece.

Warranty

Your Harvey Stethoscope comes with a five-year warranty on the binaural and chestpiece assembly. All other parts are warrantied for a period of two years against parts or manufacturing defects. During the warranty period, Welch Allyn will repair or replace the product at no charge.

Your Harvey DLX Stethoscope comes with a 10-year warranty against parts and manufacturing defects. Should a material or workmanship defect be discovered, Welch Allyn will repair or replace the product at no charge.

This warranty does not cover damage in shipment or failure caused by tampering, obvious carelessness or abuse.

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Welch Allyn Ltd. Singapore
Welch Allyn Ltd. Thailand
Medical Division, Welch Allyn Ltd.