

Excel Care[®] Bariatric Bed

Adaptable bariatric care.



Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Offering flexibility to help deliver exceptional care for a range of bariatric patient needs.

From lower acuity care to the ICU, hospitals are adjusting to accommodate an increasing population of bariatric patients.

Because increased lengths of stay, injuries related to mobility issues and higher readmission rates can be costly, it's critical to have access to the right equipment to provide appropriate bariatric care for a range of conditions.^{1,2}

The Excel Care® Bariatric Bed helps you deliver exceptional care with:

- **Flexible patient support**
- **Simplified care**
- **Support patient independence**

Hill-Rom can help **transform bariatric care.**





Flexible patient support.

The Excel Care® Bariatric Bed provides flexible support to help you manage safety issues and patient outcomes.

- Pressure redistribution surface helps treat and prevent pressure ulcers.
- Continuous air flow keeps patient skin cool and dry.
- Turn Assist helps make it easier to perform linen changes and wound checks.
- Air surface remains inflated during transport.
- Continuous Lateral Rotation Therapy (CLRT) supports protocols for helping prevent ventilator associated complications.
- Low bed height and seat deflation help make transition to and from bed easy and safe.





Simplified care.

From repositioning patients to performing other routine care, the Excel Care® Bed helps simplify the way patients are handled with easy-to-use integrated features.



- Expandable width frame and surface adjust from 40" to 50" for patient comfort.
- Adjustable bed length helps minimize patient migration and the need for patient repositioning.
- Automated features such as an integrated scale to help offer multiple options to manage your patient — in the comfort of their own bed.

Support patient independence.

The Excel Care® Bariatric Bed helps safeguard patient dignity with accessible controls that help support patient independence.



- With multiple grip points and single post patient-helper trapeze, patients can assist caregiver during repositioning or make themselves more comfortable.
- Easy-to-reach pendant offers easy access and a convenient way for patients to comfortably adjust themselves in bed.
- 1,000 lb patient weight capacity

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

1. Rice S. Retrofitting hospitals for obese patients. Modern Healthcare. 2014 Feb. modernhealthcare.com/article/20140208/MAGAZINE/302089980.
2. McGinley L. Best Practices for Safe Handling of the Morbidly Obese Patient. Bariatric Nursing and Surgical Patient Care. Volume 3, Number 4, 2008.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
134695 rev 7 02-FEB-2016 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom