



Bed Height

Low position – floor to center of deck	13" (33 cm)*
Low position – floor to edge of deck	17" (43.2 cm)*
High position – floor to center of deck	30" (76.2 cm)**
High position – floor to edge of deck	34" (86.4 cm)**
Maintains low height when placed in chair	Yes
Headboard/push handles maintain same height regardless of patient surface height	Yes

Sleep Surface

Width	35.5" (90.2 cm)
Length, fully extended	86" (218.0 cm)
Length, fully retracted	75" (191.0 cm)
Radiolucent sleep deck	No
Surface weight limit	500 lbs

Overall Width

Siderails stowed	37" (94.0 cm)
Siderails up	40" (102.0 cm)

Overall Length

Fully extended	94.5" (240.0 cm)
Fully retracted	82.5" (210.0 cm)

Bed Angles

Head section	0° to 65°
Knee section	0° to 16°
Foot section	0° to -27°
Max Trend/Reverse Trend	15° / 10°
Head and Trend angle indicators	Yes

Caster and Braking System

Optional IntelliDrive® power transport	Yes
Caster diameter	5" (12.7 cm)
with IntelliDrive® power transport	6" (15.2 cm)
Four-wheel braking	Yes
Brakes lock both rolling and swiveling	Yes

Optional Light Alerts

Feature indicators	Low low bed/chair, siderails down, brake not set, bed exit not armed
--------------------	--

Optional Bed Exit Alarm System

Sensitivity levels	3
Modes	Positioning, Exiting, Out-of-Bed, Alarm Silence

Optional Scale System

Accuracy	1% of patient weight
Weigh in any position (except Trend/Rev Trend)	Yes
Scale weight capacity	500 lbs (227 kg)

Patient Weight Limit

Patient weight, all surfaces	500 lbs (227 kg)
Safe working load, including all accessories	550 lbs (250 kg)

Patient Controls

Siderail patient bed controls	Standard
Docking pendant	Optional
SideCom® patient communication controls (Nurse call, universal TV/radio, and lighting)	Optional

Electrical Characteristics

Backup battery power	Standard
Standards test to	UL2601 & IEC 60601-1

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER 178171 rev 3 24-APR-2014 ENG – US

Zinc Number US-AC-1213-0641

www.hill-rom.com

USA 800-445-3730
CANADA 800-267-2337

*Optional IntelliDrive® power transport adds 4" (10.2 cm) height to the bed in its low position

**Optional IntelliDrive® power transport adds 1.5" (3.8 cm) height to the bed in its high position

Enhancing outcomes for
patients and their caregivers:

Hill-Rom