

# Tri-Flex II™ Bariatric Bed

Dependable Bariatric Care.



Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**

# Reliable and dependable choice for the care of bariatric patients.

Hospitals around the country face many challenges, including cutbacks, staffing shortages, and limited resources. Because of the unique medical demands, bariatric patients often require extra care, which can place even more strain on resources and staff.

The Tri-Flex II bariatric bed delivers:

- Dependable care
- Reliability
- Patient comfort



Hill-Rom can help  
**transform bariatric care.**

## Dependable Care.

- Expandable frame (37" - 48")\* for patient comfort.
- Integrated scale.
- Battery backup eliminates need to find outlet when repositioning or transferring patients.
- Easy to use caregiver controls.



## Reliability.

- Impulse Drive System\*\* allows caregiver to easily transport patient in bed, maneuvering in and out of doorways, elevators, and turns with ease.
- Low air loss and alternating pressure surface (Synergy® Air Elite) helps prevent and treat pressure ulcers.
- Low deck height of 14" and drop down rails help make the transition to and from bed easier.

## Patient Comfort.

- Patient-helper trapeze helps patients assist caregiver during repositioning or to help make themselves more comfortable.
- 1,000 lb patient weight capacity.



\* Rental (Capital: 37" – 48")

\*\*Available for capital purchase only

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
186703 rev 2 05-FEB-2016 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

**USA 800-445-3730**  
**Canada 800-267-2337**

**[www.hill-rom.com](http://www.hill-rom.com)**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**