

# Hill-Rom® Stretcher Accessories

Go Beyond Mobility



Enhancing outcomes for patients and their caregivers:



# Stretcher Accessories



## Patient Tray

product no. **P490B**

(Works on all widths of Hill-Rom® Stretchers)

Available for use on all widths of Hill-Rom stretchers without adjustment. This tray can also be placed on any location on the siderails to an area that is most accessible by the patient.



## One Button Adjustable IV Pole

product no. **P2217A**

Simply press the blue button and adjust the IV pole to desired height.



## Siderail Pads

product no. **P931BT**



## Scale Upgrade Kit

product no. **P746A01** (26")

product no. **P746A02** (30")



### Head/Footboard

product no. **P4120T**

### Chart holder available

product no. **P361**



This item can be used on any width Hill-Rom stretcher as a head or foot board. This item can also be easily removed in case of emergency, when you would need quick patient access.



### Foot Extender Pad

product no. **P929G1** (3")

product no. **P929G2** (4")  
(Requires P350T)



### Convertible Footboard

product no. **P350CTM07**

(Footboard, monitor shelf, extender)

This versatile footboard can be used as a footboard, shelf or bed extender.



**Ankle Stirrups**  
product no. **P347AT**  
(Folds away when not in use)



**Paper Roll Dispenser**  
product no. **P364AT01**  
(Standard Width 26")  
product no. **P364AT02**  
(Wide Width 30")



**Vertical Oxygen Tank Holder**  
product no. **P27601**  
(Mounts in IV sockets)



**Utility Tray**  
product no. **P297B01**  
(Standard Width 26")  
product no. **P297B02**  
(Wide Width 30")



**Liquid Oxygen Tank Holder**  
product no. **P273**




**Transport Strap**  
product no. **P349**



**Armboard**  
product no. **P344CT**



**Upright Chest X-Ray Cassette Holder**  
product no. **P279AT**



**ISS Transfer Pole**  
product no. **P158**  
(Socket Adapter P163 Required)

**Socket Adapter**  
product no. **P163**



**IV Transporter Device**  
product no. **P491**



**Lateral X-Ray Cassette Holder**  
product no. **P264**

# Surgical Stretcher Accessories



**Superior Wrist Rest**  
product no. **P262A01**




**Temporal Wrist Rest**  
product no. **P262A02**  
(Mounts on left or right side)



**Gas Delivery/  
Drape Support Device**  
product no. **P263**



**Integrated PACU Extenders**  
product no. **P261ECP**



**Head Restraint Strap**  
product no. **P449**

# OB/GYN Stretcher Accessories



**Detachable Catch/  
Fluid Basin**  
product no. **P265**



**Telescoping Calf Supports**  
product no. **P35745AT**

## Stretcher Accessory Compatibility

Accessory	Transport	Procedural	Trauma	Surgical	Electric	OB/ GYN
<b>Ankle Stirrups: P347AT</b>	X	X	X	X	X	
<b>Armboard: P344CT</b>	X	X	X		X	X
<b>Chart Holder: P361</b> (Attaches to Head/Footboard)	X	X	X	X	X	X
<b>Convertible Footboard: P350T</b>	X	X	X	X (footboard only)	X	X (headboard only)
<b>Detachable Catch/ Fluid Basin: P265</b>						X
<b>Foot Extender Pad: P929G1 (3"), P929G2 (4")</b> (Requires Convertible Footboard P350T)	X	X	X	X	X	X
<b>Gas Delivery Drape Support Device: P263</b>				X		
<b>Head/ Footboard: P4120T</b>	X	X	X	X (footboard only)	X	X (headboard only)
<b>Head Restraint Strap: P449</b>				X		
<b>Integrated PACU Extenders: P261ECP</b>				X		
<b>ISS Transfer Pole: P158</b> (Socket Adapter P163 Required)	X	X	X	X	X	X
<b>IV Transporter Device: P491</b>	X	X	X	X	X	X
<b>Lateral X-Ray Cassette Holder: P264</b> (Allows the lateral cassette holder to be mounted from head to foot)			X			
<b>Liquid Oxygen Tank Holder: P273</b>	X	X	X	X	X	X
<b>One Button IV Pole: P2217A</b>	X	X	X	X	X	X
<b>Paper Roll Dispenser: P364AT01 (26"), P364AT02 (30")</b> (Mounting point depends on configuration)	X	X	X			
<b>Patient Tray: P490B</b>	X	X	X	X	X	X
<b>Scale Upgrade Kit: P746A01 (26"), P746A02 (30")</b>		X	X		X	
<b>Siderail Pads: P931BT</b>	X	X	X	X	X	X
<b>Stretcher Pillow: P1425C</b>	X	X	X		X	X
<b>Superior Wrist Rest: P262A01</b>				X		
<b>Telescoping Calf Supports: P35745AT</b>						X
<b>Temporal Wrist Rest: P262A02</b> (Mounts onto left or right side)				X		
<b>Transport Strap: P349</b> (1 Per package)	X	X	X	X	X	X
<b>Upright Chest X-ray Cassette Holder: P279AT</b>	X	X	X	X		X
<b>Utility Tray: P297B01 &amp; 02</b>		X				
<b>Verticle Oxygen Tank Holder: P27601</b> (Mounts into IV Sockets)	X	X	X	X	X	X

# Hill-Rom® Stretcher Portfolio



HILL-ROM®  
PROCEDURAL STRETCHER



HILL-ROM®  
PROCEDURAL STRETCHER  
WITH INTELLIDRIVE®  
POWERED TRANSPORT



HILL-ROM®  
PROCEDURAL STRETCHER  
WITH ENHANCED SAFETY  
SIDERAILS



HILL-ROM®  
TRANSPORT STRETCHER



HILL-ROM®  
ELECTRIC STRETCHER



HILL-ROM®  
TRAUMA STRETCHER



HILL-ROM®  
OB-GYN STRETCHER



HILL-ROM®  
SURGICAL STRETCHER

***700 lbs — Weight capacity of all Hill-Rom® Stretchers***

**WARRANTY:**

**FRAMES & WELDS FOR LIFE OF PRODUCT**

**3 YEAR PARTS**

**2 YEAR MATTRESS**

**1 YEAR SERVICE**

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
198358 rev 2 10-FEB-2017 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

**USA 800-445-3730**  
**Canada 800-267-2337**

**[www.hill-rom.com](http://www.hill-rom.com)**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**