

# Centuris™ Pro bed with split siderails

A cost-effective solution for quality care



## Features and benefits

**AutoContour™ and SlideGuard™** minimise patient migration in bed, reducing abdominal compression during bed articulation

Easy-to-read **Line-Of-Site™ backrest indicator** assists in implementing caregiver protocols

The lowering of the **split siderails** can be activated using one hand and feature a smooth, rapid self-lowering mechanism allowing for quick access to patients

Stable support is provided during bed exit via an **ergonomically-placed side egress grip handle**

**Functional lock-out on controls** offers enhanced operational safety

A **220 kg Safe Working Load** accommodates a greater range of patients

**150 mm single band castors including steering castor** allow for ease of manoeuvrability

**Centralised braking** enhances patient safety

**The Trendelenburg/Reverse Trendelenburg function** facilitates easier patient repositioning

Choose from either a hand-held **pendant control or integrated controls**, both of which simplify operational functioning of the bed for caregivers

## Centuris Pro bed with split siderails

| Technical Specifications       |                                     |              |
|--------------------------------|-------------------------------------|--------------|
| Safe Working Load              |                                     | 220 kg       |
| Bed weight                     |                                     | 133 kg       |
| Height adjustment              |                                     | 39 to 77 cm  |
| SlideGuard sliding backrest    |                                     | up to 12 cm  |
| Protection class               |                                     | IPX4         |
| Bed angles                     | Head section                        | 0 to 65°     |
|                                | Thigh section                       | 0 to 28°     |
|                                | Foot section                        | -22 to -3°   |
|                                | Trendelenburg/Reverse Trendelenburg | -17 to +17°  |
| Sleep deck                     | Width                               | 90 cm        |
|                                | Length                              | 200 cm       |
| Siderail height                |                                     | 39.3 cm      |
| Bilateral Mechanical CPR lever |                                     | ●            |
| Steering castor                | At head end                         | ●            |
|                                | At foot end                         | ○            |
| Castors                        | 150 mm single band                  | ●            |
| Controls                       | Integrated in siderails or pendant  | ●            |
| Central brake                  |                                     | at foot end  |
| Lifting/IV pole holder         |                                     | at 4 corners |
| Bed extension and linen holder |                                     | ○            |

● = standard    ○ = optional



Integrated controls



Line-Of-Site backrest indicator



Single-hand siderail release



Pendant control

Hill-Rom is a leading global medical technology company whose products, services and more than 10,000 employees worldwide help people get better care inside and outside the hospital. Our innovations in five core areas – **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health** – improve clinical and economic outcomes and ensure caregivers in more than 100 countries have the products they need to protect their patients, speed up recoveries and manage conditions. **Every day, around the world, we enhance outcomes for patients and their caregivers.** Learn more at [hill-rom.com](http://hill-rom.com).

Not all products/options are available in all countries.  
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

The Centuris Pro bed is intended for use in intensive, general and ambulatory care or care during long hospital stays for patients with a weight  $\geq 40$  kg, height  $\geq 146$  cm and BMI  $\geq 17$ .

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. Doc. No: SEN133303-04, 6 June 2018

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**