

Approved PDAs for use with Welch Allyn[®] AcuityLink[®] Clinician Notifier (SW 1.5x)

This list indicates the Welch Allyn approved PDA devices. If your device is not on this list, contact your Welch Allyn account representative.

Installation of Clinician Notifier software on a mobile device not specified by Welch Allyn can interfere with effective Clinician Notifier program performance, including alarm notification. Only use devices approved by Welch Allyn.

Approved PDAs

For PDAs listed here:

- **Motorola ES400**
- **Motorola MC75A**

Use the following settings:

- Operating system: Windows Mobile 6.5x (see device documentation)
- Screen resolution: 640 x 480
- For information on installing Clinician Notifier on these devices, see the instructions under the Technical tab at <http://www.welchallyn.com/apps/products/product.jsp?id=16-vo-96-1232546273407> and click on the device name.