

Approved PDAs for use with Welch Allyn® AcuityLink® Clinician Notifier (SW 1.1x-1.4x)

This list indicates the Welch Allyn approved PDA devices. If your device is not on this list, contact your Welch Allyn account representative.

Installation of Clinician Notifier software on a mobile device not specified by Welch Allyn can interfere with effective Clinician Notifier program performance, including alarm notification. Only use devices approved by Welch Allyn.

Approved PDAs

Motorola MC55/Symbol MC5590

- Operating system: Windows Mobile 6.1 classic
- Screen resolution: 320 x 240

Motorola MC50/MC5040 (upgrade only)

- Operating system: Windows mobile 2003 SE 4.21.1088
- Screen resolution: 320 x 240
- For information on upgrading Clinician Notifier on this PDA, see the instructions under the Technical tab at <http://www.welchallyn.com/apps/products/product.jsp?id=16-vo-96-1232546273407>, and click on the device name.

For PDAs listed here:

- **HHP 7900**
- **HHP 9500**
- **Honeywell Dolphin 7600**
- **Honeywell Dolphin 9900**
- **Motorola MC7090**
- **HP iPaq 111**
- **Pharos 600e**
- **Unitech PA600 MCA**

Use the following settings:

- Operating system: Pocket PC 2003 SE, Windows Mobile 5.x, Windows Mobile 6.0x or 6.1x (see device documentation)
- Screen resolution: 320 x 240
- For information on installing Clinician Notifier on these devices, see the instructions under the Technical tab at <http://www.welchallyn.com/apps/products/product.jsp?id=16-vo-96-1232546273407> and click on the device name.