

# Welch Allyn Scale-Tronix® Stow-A-Weigh® Scales



Scale-Tronix Stow-A-Weigh scales deliver a big impact with a small footprint. These compact scales provide the fast, dependable weight measurements you need to help improve patient care and outcomes while saving on floor space.

These innovative fold-up scales are ideal for clinics, triage areas, exam rooms or anywhere space is limited. With stand-on, wheelchair and stretcher models available, it's easy to find the right low-profile weighing solution for your environment.

Scales that are big on performance and small on space—it's easy to see why Scale-Tronix is a trusted brand for healthcare environments.



## Durable. Long Lasting. Designed for Demanding Healthcare Environments.



### 5202 Stow-A-Weigh Stand-On Scale

- Mounts snug to wall with slender 4" profile when not in use
- Weighing platform is counter-balanced to fold down easily—no bending or lifting required
- Low-profile platform helps reduce patient fall risk

### 6202 Stow-A-Weigh Wheelchair Scale

- Easily accommodates most wheelchairs
- Three scales in one: wheelchair, chair and stand-on

### 6202D Dual Platform Stow-A-Weigh Wheelchair Scale

- Easily accommodates most stretchers
- Four scales in one: stretcher, wheelchair, chair and stand-on

## Find the Right Scale for You

	5202 Stow-A-Weigh Stand-On Scale	6202 Stow-A-Weigh Wheelchair Scale	6202D Stow-A-Weigh Scale for Stretchers, Wheelchairs, Chairs, Scooters and Ambulatory Patients
<b>Weighing Capacity</b>	660 lb / 300 kg	660 lb / 300 kg	880 lb / 400 kg
<b>Dimensions</b>	Platform: 18" w x 18" d x 1.5" h Overall: 27" w x 4" d x 27" h	Platform: 24" w x 30" d x 1.5" h Overall: 40" w x 4" d x 40" h	Platform A: 24" w x 30" d x 1.5" h Overall: 40" w x 4" d x 40" h Platform B: 12" w x 30" d x 1.5" h Overall: 27" w x 4" d x 40" h
<b>Resolution</b>	0.1 lb / 100 gm	0.1 lb / 100 gm	0.5 lb / 200 gm
<b>Power Source (Cordless)</b>	6 D-size disposable alkaline batteries	6 D-size disposable alkaline batteries. Power adapter available	6 D-size disposable alkaline batteries. Power adapter available
<b>Readout</b>	Digital LED display in pounds and kilograms, or kilograms only		
<b>Optional Accessories</b>	Height gauge (wall mount, stainless steel, in/cm)	Height gauge (wall mount, stainless steel, in/cm), printer, AC power adapter, thermal printer paper (15 rolls/box)	
<b>Automatic Zero</b>	Automatically returns to zero, ready for next patient		
<b>Weight Recall</b>	Recalls last weight		
<b>Reweigh</b>	Recomputes patient's weight while patient is on scale		
<b>Computer Connectivity</b>	N/A	RS 232 output (optional)	

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361