

WelchAllyn®



Welch Allyn
MacroView™ Otoscope




A bigger, better, sharper view

Enhancing visualization of the ear, helping you perform routine exams with **ease and confidence.**



Integrated Wall System



Diagnostic Desk Set




Diagnostic Set

Better clarity and definition of landmarks

MacroView optics increase **magnification by 30%** versus a traditional otoscope, giving you a nearly full view of the tympanic membrane. Coupled with the widest focus range, MacroView offers a **bigger, better, sharper view.**

2X FIELD OF VIEW



Standard Otoscope

MacroView Otoscope


Welch Allyn LED Lamps

LED technology provides the benefits of LED while operating at a color similar to halogen.




WITH THE FAMILIARITY OF HALOGEN

SureColor™ LED



Optimized color temperature, similar to what you are used to with halogen

Most Competitors' LED



Light too bright or too white can wash out critical details

LED IS ENERGY EFFICIENT

3X the on-time vs. halogen with Welch Allyn rechargeable handles

Lower operating temperature intended for greater patient/clinician comfort



HALOGEN



LED



7-YEAR WARRANTY

Welch Allyn MacroView Otoscope



Proprietary **"Tip Grip"** helps ensure ear speculum is secured and ejects used ear speculum for easy disposal

Default focus position; ideal for most ear examinations

Adjustable focus wheel helps you accommodate for variable ear canal lengths and farsighted eyes, zooming in or out to fine tune view

Optional throat/nasal illuminator adds one more tool to your practice; easy to adjust on 3.5 V handle

3.5V Lithium Ion rechargeable handle (shown) provides 2.4X more power compared to our Nickel Cadmium rechargeable handle

ACCESSORIES



52432-U
Kleenspec® Disposable
Ear Specula



52700
Slotted Specula
Instrumentation Tip



23804
Otoscope Insufflation
Bulb and Tube



24330
Sofseal Specula Tip

Welch Allyn Digital MacroView™ Otoscope

Supports the diagnosis with better communication
to enhance patient compliance

- Store digital images in patient records to support documentation and diagnosis
- Offers the ability to view clear, live and still images of small details such as tubes and hearing aids
- An excellent training aid to support clinical education



Configurations

77791-MX	Welch Allyn Green Series™ 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope, KleenSpec Disposable Specula Dispenser, Wall Aneroid Sphygmomanometer and SureTemp® Plus 690 Electronic
77792-M	Welch Allyn Green Series 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope, KleenSpec Disposable Specula Dispenser and Wall Aneroid Sphygmomanometer
77792-MNOBP	Welch Allyn Green Series 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope and KleenSpec Disposable Specula Dispenser
71641-MS	Welch Allyn 3.5 V Halogen HPX™ Diagnostic Set including MacroView Otoscope, Coaxial Ophthalmoscope, two Lithium Ion Rechargeable Handles and Universal Desk Charger
71641-M	Welch Allyn 3.5 V Halogen HPX Diagnostic Set including MacroView Otoscope, Coaxial Ophthalmoscope, two Nickel Cadmium Rechargeable Handles and Universal Desk Charger
23920	Welch Allyn 3.5 V Halogen HPX Digital MacroView Fiber-Optic Otoscope; Mini-USB Cable (10.0 ft/3.0 m) and set of Reusable Ear Specula
Accessories	
23804	Insufflation Bulb
52700	4.00 mm KleenSpec Slotted Disposable Specula Instrumentation Tip
52432-U	2.75 mm Universal KleenSpec Pediatric Disposable Ear Specula
52434-U	4.25 mm Universal KleenSpec Disposable Ear Specula
24320	SofSeal Specula Tip for Disposable Ear Specula, small
24330	SofSeal Specula Tip for Disposable Ear Specula, medium

For more information on the MacroView Otoscope, please contact your Welch Allyn sales representative or visit www.welchallyn.com.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361



www.welchallyn.com

WelchAllyn®