



**Hillrom™**

# ADVANCING CONNECTED CARE™ WITH ECG CONNECTIVITY

## Hillrom and Cerner

We're helping you improve clinical workflows and reduce manual documentation for ECG data by sending it right to your EHR. See how the Welch Allyn® Diagnostic Cardiology Suite™ ECG can help create a workflow that works best for your practice.

## Unconnected Workflow

- 1 Order an ECG
- 2 Manually document patient demographics in the ECG software
- 3 Capture ECG and print test results for physician review
- 4 Add physician notes to ECG printout
- 5 Scan the ECG printout and import into Cerner Millennium® EHR
- 6 Add diagnosis in Cerner Millennium EHR
- 7 Complete order

## Connected Workflow


- 1 Order ECG in Cerner Millennium EHR
- 2 Launch Diagnostic Cardiology Suite ECG
- 3 Capture ECG and send to Cerner ECG Management
- 4 Physician opens ECG in Cerner ECG Management to review, edit and sign
- 5 Image of the final ECG result and discrete data are automatically sent to Cerner Millennium EHR

## Connectivity Counts

Help save time and improve data accuracy by streamlining your workflows.

- Save time previously spent manually entering data into the EHR
- Improve accuracy by automatically associating results with the EHR orders
- Work more efficiently by freeing up your time for other patient priorities

Automate your workflow with the Diagnostic Cardiology Suite ECG by connecting directly to Cerner Millennium EHR.



Contact your Cerner representative to learn more.



[hillrom.com](http://hillrom.com)

\*A Wireless Acquisition Module or an AM™12 Module (via a USB connection) can be used with the Diagnostic Cardiology Suite to capture an ECG.

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.

© 2019 Welch Allyn, Inc. ALL RIGHTS RESERVED. MC15250 80022623 Ver A 2019-09