



**Hillrom™**



# SHORT TRANSFERS MADE SIMPLE

## HANDYSHEET™ AND HANDYTUBE™ MANUAL TRANSFER AIDS.

Short transfers are often anything but simple. That's why we designed the HandySheet and HandyTube manual transfer aids to help you reposition patients and apply slings safely—and with less strain.

- Slide patients more easily with high-quality fabric that reduces friction between the patient and the underlying surface.
- Reposition a wide range of patients - from normal weight to patients of size - as well as support prone positioning.
- Expect long-term performance with fabric that holds up to repeated washing.



#### HANDYSHEET™

The HandySheet is an excellent aid for applying a sling, especially on patients of size, pain-sensitive patients and patients with contractures. Designed for use in pairs, it comes in a two-sheet set and is available in two sizes:

- Long versions have bed attachment loops so fewer caregivers are needed.
- Short versions have small pockets to help you apply a sling in small spaces (e.g., narrow wheelchairs).



#### HANDYTUBE™

The HandyTube is a sliding sheet designed to help you move patients higher up in bed, turn them in bed or transfer them between beds.

It is available in short, long and wide models to meet a variety of repositioning needs.

### ORDERING INFORMATION

Part Number	Description	Size	Quantity
3715015-10	HandySheet Short, without handles, size M	70 x 70 cm (28 x 28 in)	Sold in sets of 10 pairs
3715017-10	HandySheet Short, without handles, size XL	120 x 85 cm (47 x 34 in)	Sold in sets of 10 pairs
3716015-10	HandySheet Long, with handles, size M	90 x 185 cm (35 x 75 in)	Sold in sets of 10 pairs
3716017-10	HandySheet Long, with handles, size XL	140 x 215 cm (55 x 85 in)	Sold in sets of 10 pairs
3710015-10	HandyTube Short, without handles, Size M	70 x 90 cm (28 x 35 in)	Sold in sets of 10 pcs
3711015-10	HandyTube Long, without handles, Size M	90 x 185 cm (35 x 73 in)	Sold in sets of 10 pcs
3712015-10	HandyTube Wide, with handles, Size M	140x140 cm (55 x 55 in)	Sold in sets of 10 pcs

For more information, please contact your local distributor or Hillrom sales representative.

**hillrom.com**

The HandyTube™ and the HandySheet™ are developed by Hill-Rom to facilitate patient transfers while in a prone position. They are manufactured from a low-friction material that minimizes friction between the patient and the bedding. HandyTube and HandySheet are not intended for lifting or transferring patients from a lying to sitting position.

Class I

Manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå, Sweden

This device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.