



ACCELLA™ SMART BED

Designed for patients and caregivers specifically to solve challenges across mid-to-high acuity care settings to help improve patient safety, enhance outcomes, and optimize workflow.

Improve Patient Safety

- 3-mode Bed Exit Alert*
- Intelligent Nightlight
- Bed in Lowest Position Indicator
- Audible Brake Off Indicator
- Steam Clean Validated

Enhance Outcomes

- Patient Migration Management with SlideGuard™ Technology
- Intelligent AutoContour™
- Pressure Redistribution
- Accella Therapy™ Integrated Surface w/Advanced Microclimate® Technology*
- One-Touch Easy Chair™

Optimize Workflow

- Simple Weigh/BMI In-Bed Scale
- Graphical Caregiver Interface (GCI)
- Auto CPR/Lever
- Head-Of-Bed (HOB) Angle Alert
- One-Touch Side Egress
- Wireless Smart Connectivity Module

* Optional Accella™ Smart Bed features and enhancements.



Hillrom™

Versatile and innovative, the Accella™ Smart bed offers flexible solutions across a range of mid-to-high acuity care settings and provides a future-proof design for future advancements.



TECHNICAL SPECIFICATIONS

Safe working load		250kg
Bed weight		170kg
Height adjustment		38.6 cm to 80.8 cm
SlideGuard™ technology sliding backrest		up to 12 cm
Protection class		IPX4
Bed angles	Head section	0 to 65°
	Thigh section	0 to 28°
	Foot section	-22° to -3°
	Trendelenburg/ Reverse Trendelenburg	-17 to +17°
Sleep deck	Width	90 cm
	Length	200 cm
Sliderail height		39.3 cm

Optional Accella™ Smart Bed features and enhancements include:

- 3-mode Bed Exit Alert
- Accella™ Therapy Integrated Surface
- Integrated X-Ray
- Bed Extension Linen Holder
- Patient Pendant
- Double Band Castor Options – (125mm/150mm)
- 5th Castor
- Brake/Steer Pedal at Head and Foot
- Nurse Call

For further information about our products and services, please contact your local Hillrom representative or visit our website.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

The Hill-Rom 900 Accella bed is intended for use in intensive, acute and ambulatory care settings for patients with a weight ≥ 40 kg, height ≥ 146 cm and BMI ≥ 17 . Medical device (93/42/EEC): Class I, Weigh system (2014/31/EEC): Class III, Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France
This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.