



Hillrom™



The Hillrom TS 7000dV Surgical Table is designed to work with Integrated Table Motion for the *da Vinci Xi*® Surgical System from Intuitive Surgical.

TS 7000DV SURGICAL TABLE

REDEFINING PATIENT POSITIONING IN DA VINCI® SURGERY

Integrated Table Motion (ITM) enables the *da Vinci Xi*® System to connect to Hillrom's TS 7000dV Surgical Table so the patient can be dynamically positioned while the surgeon operates.

Iso-centered motion is achieved through ultra-smooth, electromechanical actuation allowing for simultaneous Trendelenburg or reverse Trendelenburg and slide creating a virtual pivot point.



TS 7000dV Options

The TS 7000dV Surgical Table is outfitted with patient positioning accessories for a complete solution. Package options may vary, though mainly consist of the following accessories:

- Double Joint Head Section
- Single Part Leg Section
- dV Remote Control
- dV Connection Cable
- Lightweight Transfer Board
- Armboards
- Component Cart
- Patient Restraint Strap
- Rotatable Footboard & Clamps
- Multi-task Armboard & Clamp*

With any TS 7000dV Operating Table two years of comprehensive service coverage and maintenance is included.



Dramatic Range of Motion

Patient Weight Capacity	1,000 lb (450 kg)	Exceptional weight capacity in any articulation, including slide
Height Range ¹	23 - 45" (59 - 114 cm)	Ultra low and high platform for the surgical team, including the surgical robot
Trendelenburg / reverse Trendelenburg	± 45°	Superior positioning to expose target anatomy, including during robotic-assisted surgery

Additional Accessories

Hillrom offers an additional accessory package of specialty components that optimize patient positioning capabilities for complex robotic-assisted surgeries.



Two-Part Leg Section for midline access.



Pelvic Extension for enhanced imaging capability.



TrenGuard for secure, steep Trendelenburg positioning.

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

*Not included in all packages; Packages subject to change

da Vinci and da Vinci Xi are trademarks or registered trademarks of Intuitive Surgical, Inc.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR76501 rev 1 16-APR-2020 ENG – US