



Hillrom™



TS 7500 OR Table with floor mounted column and U24 tabletop;
TS 7500 OR Table wireless remote (below).

UNIVERSAL TABLETOPS

SPECIALTY SURGICAL TABLETOPS FOR A WIDE RANGE OF PROCEDURES

The TS 7500 Operating Table offers a high number of positioning options for a wide variety of applications: from general and neurosurgery to gynecology, urology, traumatology, and orthopedics. Flexibility is further increased by an extensive range of positioning pads and special table accessories.





BARIATRIC



GENERAL SURGERY



NEUROSURGERY



SPINAL SURGERY



THORACIC SURGERY



ORTHOPEDICS



GYNECOLOGY



ORTHOPEDICS AND TRAUMA



TRANSFER TABLETOPS



TECHNICAL DETAILS

	U14 Tabletop	U24 Tabletop	U26 Tabletop
Material no.	# 1846376	# 1846381	# 1846386
Dimensions (L x W)	32.8" x 23.6" (832 mm x 600 mm)	42.1" x 23.6" (1070 mm x 600 mm)	42.1" x 23.6" (1070 mm x 600 mm)
Height range	20.5" – 42.1" (520 mm – 1070 mm)	20.5" – 42.1" (520 mm – 1070 mm)	20.5" – 42.1" (520 mm – 1070 mm)
X-ray area	14.7" (370 mm)	Head end: 27.0" (685 mm) Foot end: 14.7" (370 mm)	Head end: 27.0" (685 mm) Foot end: 14.7" (370 mm)
Longitudinal shift	18.1" (460 mm)	15.7" (400 mm)	15.7" (400 mm)
Trendelenburg / reverse Trendelenburg	± 45°	± 45°	± 45°
Tilt	± 30°	± 30°	± 30°
Maximum patient weight	880 lbs (400 kg)	880 lbs (400 kg)	880 lbs (400 kg)

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.