

# PROGRESSA<sup>®</sup> SMART BED SYSTEM

Empower patients and caregivers with in-bed therapies to help get your most critical patients better sooner.



## STANDARD FEATURES TO HELP OPTIMIZE PATIENT SAFETY AND IMPROVE CAREGIVER SATISFACTION

|  |
|--|
| SlideGuard™ Frame Technology                     |
| Patient Controls Backlighting                    |
| FlexAfoot™ Retractable Foot Control              |
| Advanced MicroClimate® Technology                |
| 30°/45° Head of Bed Alarm                        |
| In-Bed Scale                                     |
| Line Managers                                    |
| Three Mode Bed Exit Alarm                        |
| Caregiver Pendant                                |
| Caregiver Foot Controls                          |
| Night Light                                      |
| Integrated Oxygen Tank Holder                    |
| Accessory Outlet                                 |
| Hands Free Emergency CPR & Trendelenburg         |
| NaviCare Ready - Local Alarm On/Off              |
| Point of Care™ Siderails Control                 |
| Both Siderails with Point-of-Care Touch Screen   |
| Smart Bed Ready (Sidecom Required)               |
| Obstacle Detect™ System at Both Sides of the Bed |
| Drainage Bag Holders                             |
| Four IV Sockets                                  |

## OPTIONAL FEATURES TO SUPPORT THERAPY-FOCUSED TECHNOLOGIES

|  |
|--|
| Prevention Surface                               |
| Therapy Surface                                  |
| Pulmonary Surface                                |
| Continuous Lateral Rotation Therapy (CLRT)       |
| Percussion & Vibration Therapy (P&V)             |
| Nurse Call, Univ TV & Lighting                   |
| Siderail Nurse Call                              |
| Bed Exit Alarm to Nurse Call                     |
| 5th Wheel Steering                               |
| IntelliDrive® Powered Transport Steering         |
| Transport Shelf                                  |
| Wireless Connectivity                            |
| Permanent IV Pole                                |
| Full Chair Egress                                |
| StayInPlace™ Patient Migration Management System |
| WatchCare™ Incontinence Management System        |



Hillrom™

# Progressa® Smart Bed System

## TECHNICAL SPECIFICATIONS

|  |  |
|--|--|
| Total Length                             |  |
| Foot Section Extended                    | 98" (248.9 cm)                                     |
| Foot Section Retracted                   | 88" (223.5 cm)                                     |
| Maximum Width                            |  |
| Side Rails Down                          | 40.5" (103 cm)                                     |
| Side Rails Up                            | 40.5" (103 cm)                                     |
| Maximum Side Rail Height                 |  |
| Above Sleep Deck                         | 17" (43.2 cm)                                      |
| Highest Position                         |  |
| Top of Seat Section to Floor             | 35.7" (90.7 cm)                                    |
| Lowest Position                          |  |
| Top of Seat Section to Floor             | 16.5" (41.9 cm)                                    |
| Minimum Underbed Clearance               | 4.3" (10.9 cm) w/o IntelliDrive® Powered Transport |
| Standard Caster Size                     | 6" (15.2 cm)                                       |
| <b>Articulations</b>                     |  |
| Head Section                             | 0° - 67°<br>77° with Chair Egress                  |
| Thigh Section                            | 0° - 30°   |
| Foot Section                             | 0° - 45°   |
| Trendelenburg Position (maximum)         | 75° with Chair Egress                              |
| Reverse Trendelenburg Position (maximum) | 13°<br>18°   |

|                                |  |
|--------------------------------|--|
| Safe Working Load              | 650 lb (295 kg)                        |
| Maximum Patient Weight         | 500 lb (227 kg)                        |
| C-Arm Clearance                | 7 - 15" (17.8-38.1 cm)                 |
| Cleanability/Infection Control | IPX4                                   |
| Radiolucent Sleep Deck         | Standard 17.7" x 23" (45 cm x 58.5 cm) |
| Voltages                       | 100/110/115/120/127/220/230/240        |
| <b>Mattress</b>                |  |
| <b>Prevention Surface</b>      |  |
| Mattress Width                 | 35" (88.9 cm)                          |
| Mattress Length                | 84" (213 cm)                           |
| Mattress Height                | 7.1" (18.1 cm)                         |
| Mattress Weight                | 31 lb (14.1 kg)                        |
| <b>Therapy Surface</b>         |  |
| Mattress Width                 | 35.5" (90.2 cm)                        |
| Mattress Length                | 84" (213 cm)                           |
| Mattress Height                | 8" (20.32 cm)                          |
| Mattress Weight                | 45 lb (20.5 kg)                        |
| <b>Pulmonary Surface</b>       |  |
| Mattress Width                 | 35.5" (90.2 cm)                        |
| Mattress Length                | 84" (213 cm)                           |
| Mattress Height                | 8" (20.32 cm)                          |
| Mattress Weight                | 48 lb (21.8 kg)                        |

## SURFACE & FRAME CONFIGURATION OPTIONS

| PROGRESSA® PREVENTION  | PROGRESSA® THERAPY   | PROGRESSA® PULMONARY   |
|--|--|--|
| Prevention Surface   | Therapy Surface  | <b>Choose one:</b><br>- Pulmonary Surface (CLRT only)<br>- Pulmonary Surface (CLRT & P&V)  |
| <b>Choose one:</b><br>- Dining Chair® Position<br>- FullChair® Position with Chair Egress<br>- FullChair® Position with Chair Egress and StayInPlace™ Technology | <b>Choose one:</b><br>- Dining Chair® Position<br>- FullChair® Position with Chair Egress<br>- FullChair® Position with Chair Egress and StayInPlace™ Technology | <b>Choose one:</b><br>- Dining Chair® Position<br>- FullChair® Position with Chair Egress<br>- FullChair® Position with Chair Egress and StayInPlace™ Technology |



Hillrom™

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR79203 rev 1 01-JUN-2020 ENG - US