



Introducing the Welch Allyn Single-Use Sheath

for the Welch Allyn KleenSpec®
Cordless Illuminator.



Now you can help prevent cross-contamination and keep your cordless illuminator even cleaner with the single-use sheath designed specifically for the Welch Allyn KleenSpec Cordless Illuminator used with the 590 series Disposable Vaginal Specula.

Simply pull a disposable sheath from the box, place it over the top of the illuminator and pull down. Then insert the illuminator into the disposable 590 series vaginal speculum and discard the sheath after your exam.



Ordering Information

Part # 59010 Disposable Sheath (100/box, 5 boxes/case)

For more information, contact your local
Welch Allyn representative or call 800.535.6663

WelchAllyn®

Advancing Frontline Care™