

Envision[®]


E700 Therapy Surface



The Envision surface represents one of the most advanced wound therapies Hill-Rom offers in a powered air surface. Differentiated solutions promote improved outcomes for patients with moderate to severe pressure injury risk:




PressureSmart[®]
Sensing Technology



Advanced Microclimate[®]
Technology




ShearRelief[®]
System



Vertical Air Cell
Technology

Envision® E700 Therapy Surface



Technical Specifications

Mattress Width (inflated)	35.5" ± 0.5" (90.2 cm ± 12.7 mm)
Mattress Length (inflated)	84.5" ± 0.5" (214.6 cm ± 12.7 mm)
Mattress Height (inflated at maximum pressures for a step deck bed such as the VersaCare® Bed or TotalCare® Bed System)	10.875" ± 1" (27.623 cm ± 3 cm)
Mattress Height (inflated at maximum pressures for a flat deck bed such as the Advanta™ Bed or Advance® Series Bed)	8" ± 1" (20 cm ± 3 cm)
Air Supply Unit Dimensions	14.5" W x 13.75" H x 8.5" D
Air Supply Unit Weight	21 lbs (10 kg)

Power Cord Length	168" +/- 1" (427 cm +/- 3 cm) See user manual
Power Requirement	120 V AC, 60 Hz, 3A
Recommended Therapeutic Patient Weight Range	70 lbs to 400 lbs (32 kg to 181 kg)

Order Numbers

Envision® E700 Wound Surface Mattress Replacement System	P741A
Envision® E700 Wound Surface on VersaCare® Bed	P741A01

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
207172 rev 1 18-DEC-2017 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom