

Mortara Web Upload

SAFE AND SECURE HOLTER AND ABP DATA TRANSFER

PRODUCT FEATURES

- **Unlimited Number of Client Sites**

Mortara Web Upload securely and automatically installs on client computers. A web-based application means no limitation to the number of remote client sites. Plus, client updates are automatic and transparent.

- **Data Security**

Secure transfer (SSL) with authentication to protect patient's privacy facilitates regulatory and HIPAA compliance. In the event of transfer interruptions, a notification is displayed and remaining data can be uploaded at a later time.

- **Simple Operation**

Wizard-like user interface allows for quick and easy H3+™, H12+™, and ambulatory blood pressure recorder data transfer, patient demographic entry, diary event annotation, and recorder media preparation. Demographic data fields may be auto-populated through an orders interface.

- **Rapid Results**

Uploaded recordings can be displayed and filtered as recordings are imported at the Hscribe™ Holter and Ambulo™ 2400 HDS analysis systems.

- **Flexible Customization**


Web server administrators can configure their own banner logo along with languages, server settings, server lists, and data storage locations.

- **Simple Installation**

Automatic web server software installation, plus web link for easy client installation and upgrades, minimizes setup and maintenance.



 **Mortara**



**MORTARA INSTRUMENT, INC.**

7865 North 86th Street
Milwaukee, WI 53224
U.S.A.

Tel: 414.354.1600
Tel: 800.231.7437
Service: 888.MORTARA
Fax: 414.354.4760

MORTARA INSTRUMENT EUROPE S.R.L.

(European Headquarters)
Via Cimarosa 103/105
40033 Casalecchio di Reno (BO)
Italy
Tel: +39.051.298.7811
Fax: +39.051.613.3582

MORTARA INSTRUMENT AUSTRALIA

PO Box 7568
Baulkham Hills NSW 2153
Unit 28, 9 Hoyle Avenue
Castle Hill NSW 2154
Australia
Tel: +61 2 8070 9303
Fax: +61 2 9899 9478

MORTARA INSTRUMENT UK LTD

Units 11 & 12, Scion House
Stirling University Innovation Park
Stirling FK9 4NF
Scotland
Tel: +44.1786.444980
Fax: +44.1786.446630

www.mortara.com

ISO 13485 CERTIFIED

**WARRANTY + SERVICE**

Mortara Instrument is committed to the highest level of customer support. Please contact us for the program which best suits your needs.

MORTARA™ WEB UPLOAD

Mortara Web Upload software is composed of two components:

- 1. Web Services** handles requests received from the client upload application.
- 2. Client Web Upload** allows clients to install the upload application on their computer from the host site web page. Clients need to visit this web page only once.

WEB SERVICES	REQUIREMENT*
WEB SERVER:	
OPERATING SYSTEM	Microsoft® Windows® 2008 or 2012 server R2 64-bit
RAM	1 GB memory (minimum)
PROCESSOR	1.6 GHZ Pentium IV (minimum)
ASP.net	Microsoft .NET Framework V2.0
INTERNET INFORMATION SERVICES	Microsoft IIS 8 or later
CLIENT WEB UPLOAD	
COMPUTER:	
OPERATING SYSTEM	Windows 7 Professional 32-bit or 64-bit and Windows 8.1 64-bit
RAM	1 GB memory (minimum)
VIDEO	1024 x 768 resolution (minimum)
ASP.net	Microsoft .NET Framework V2.0
WEB BROWSER	Microsoft Internet Explorer V8.0 or later
H12+ CF CARD UPLOAD	H12+ CF card reader (PN 9903-013)
H3+ RECORDER UPLOAD	H3+ USB interface cable (PN 25019-006-50)
AMBULO 2400 RECORDER UPLOAD	Mini USB cable (PN 008-0008-01)

*Requirements subject to change without notice.

Ambulo™, Burdick®, H3+™, H12+™, Hscribe™, Mortara® and Quinton® are trademarks or registered trademarks of Mortara Instrument, Inc. Microsoft® and Windows® are registered trademarks of Microsoft Corporation.

2015-11