

Introducing Surestream™ CO₂ Sampling Lines

Now available for improved CO₂ sampling



Oridion is pleased to offer a new line of CO₂ sampling products. Surestream sampling lines are a combination of Oridion and Welch Allyn monitoring technologies for use with intubated and non-intubated patients.

Surestream incorporates a unique sample technology for enhanced waveform integrity and accuracy that works with Welch Allyn Propaq® and Atlas™ monitors:

- **Nafion® for high humidity environments**
- **Includes both a hydrophobic fiber and pellet to prevent moisture from entering the monitor while maintaining laminar flow**

Surestream™

etCO₂ sampling for non-intubated patients

Surestream etCO₂ sampling lines for non-intubated patients have an Oridion innovative patient interface technology that reliably samples etCO₂ in both oral and nasal breathers providing an accurate, reliable and easy to use assessment of your patient's breathing.

Features and Benefits

- **Increased surface area of oral sampling prong** – Provides greater accuracy in the presence of low tidal volume
- **Uni-junction™ patented design** – Assures sampling from either the oral or nasal breathing source
- **O₂ delivery system** – Provides comfortable delivery up to 5 l/m without dilution of CO₂ sample; delivery via small pin-holes
- **Selection of patient interfaces** – Meets specific needs

Oral/Nasal Sampling with O₂ delivery



Smart SureLine™ Plus O₂
010981



Smart SureLine O₂ Pediatric
010982

Endoscopy Solution



*Smart SureLine Guardian™
012533, 012534**
shown above

Gum Comfort Pad
012541

Hook and Loop Strap
0542

Nasal Sampling with O₂



SureLine O₂ Adult
010979
shown above

SureLine O₂ Pediatric
010980

Nasal Sampling



SureLine Adult
010976
shown above

SureLine Pediatric
010977

*Pioneering CO₂ sampling system for upper endoscopy procedures. The Smart SureLine Guardian O₂ allows continuous CO₂ monitoring and O₂ delivery without compromise of patient position or endoscopic procedure.

** This product has extra long tubing (4M/13ft)

etCO₂ sampling for intubated patients

The design of Sure VentLine™ Sets effectively monitors intubated patients even in very humid environments.

Features and Benefits

- **2-port design** – Reduces clogging from moisture and secretions
- **Airway adapter** – Performs in any position
- **Available with additional Nafion** – Provides superior moisture handling
- **Able to connect to 15mm or 22 mm ET tube**
- **Effective dead space: < 5.6 cc**
- **Low dead space (7.3cc)**



Sure VentLine Set
010986



Sure VentLine H Set
(extra Nafion for humid environments)
010987

All Surestream breath-sampling products are for single patient use only and are compatible in the MRI environment.

All part numbers are for 25/pack.

Standard sampling line tubing length: 2M (6.5 ft)

Specifications are subject to change without notice.

Surestream products are latex free.

Surestream is a trademark of Oridion Medical 1987 Ltd.



Oridion[®]

Oridion Capnography Inc.
160 Gould Street
Suite 205
Needham MA 02494
USA
Tel: (888) ORIDION
Fax: (781) 453-2722
www.oridion.com

Oridion Medical 1987 Ltd.
7 Hamarpe Street
97774
Jerusalem
Israel
Tel: +972 2 589 9158
Fax: +972 2 586 6680

© 2010 Oridion Medical 1987 Ltd. All rights reserved. MM0511B