

Welch Allyn, Inc
 4341 State Street Road
 P.O. Box 220
 Skaneateles Falls, NY 13153-0220
 USA
 Telephone: 315-685-4100



December 14, 2015

To: Welch Allyn US Distributors and Customers

Re: Pro4000 and Pro6000 Probe Cover Compatibility

From: Jamie Wood, US/C Product Manager


With the introduction of the new Braun Pro6000 Ear Thermometer, we have also launched a new Probe cover for the Pro6000 device. This probe cover is also backwards compatible to the Braun Pro4000 device when proper loading technique is used (see diagram below):

The Braun Pro4000 Probe covers are only compatible with the Braun Pro4000 device and will not work on the Pro6000.

PRO 4000 Probe Cover	PRO 6000 Probe Cover
Flat Flange	Raised Flange
Compatibility: Pro 4000 <u>only</u>	Compatibility: Pro6000 and Pro4000
Part #: 05075-005	Part #: 06000-005

Follow these loading instructions for proper loading of probe covers on Braun Pro4000:


BRAUN
PRO 4000
Proper Probe
Cover Installation




1.




2.




3.



4.





This product is manufactured by Kaz USA, Inc. under license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany. ©2014 Kaz USA, Inc.